


Examination Record

General		Wound	No
Diabetes	Yes	Alcohol	No
Smoker	No	Dyslipidemia	Yes
Walking Pain	Yes	Arm pressure	120
Rest pain	Yes		


**60 year male. Diabetic.
Incompressible Arteries.
Suspected Mediasclerosis**

Pressure values and ABI/TBI




Recommended Values	
ABI	$\leq 1,40-0,90 \geq$
TBI	$\geq 0,70$
Ankle Pressure	$\geq 70\text{mmHg}$
Toe Pressure	$\geq 50\text{mmHg}$



Curves Right Ankle



Curves Left Ankle



Curves Right Toe



Curves Left Toe

